


## Complete Rugged NaI Detector Probes for In-Field Gamma Spectroscopy Measurements



- Rugged integrated NaI detectors in tough aluminum housing — available in 2" x 2" and 3" x 3" sizes.
- Integral preamplifier and high voltage (HV) for convenience and safety.
- Low power consumption (240 mW), ideal for portable applications.
- Single-cable connections for most applications.
- 2BY2-DD and 3BY3-DD available for use with ORTEC digiDART and digiDART-LF.

### Everything You Need in One Compact Package

2BY2 or 3BY3 is a complete NaI(Tl) probe, ready to go to work in your application. Each model contains a NaI(Tl) crystal (2" x 2" or 3" x 3") with photomultiplier tube (PMT) and a PMT base with integral preamplifier and HV supply securely encased in a compact aluminum housing.

### Ultra Low Power; the End of Bulky, External, HV Supplies

The model 296 ScintiPack™ PMT base includes a low-power, adjustable high-voltage supply, an active bias network, and a charge-sensing preamplifier. Incorporating the bias supply in the PMT base eliminates high voltage cable connections to bulky, external, HV supplies. The ACTIVE bias network permits high-rate operation with minimal peak shift. ScintiPack's ultra-low power consumption (240 mW) makes it attractive for portable applications.

### Single-Cable Operation

A single cable supplies power to the preamplifier and the HV supply as well as the signal out of the preamplifier. For the "-DD" versions, the connection is direct to the digiDART or digiDART-LF (also other ORTEC DSP spectrometers). An optional break-out cable can be used with the model 2BY2 or 3BY3 should you want the signal on a BNC output connector.

# 2BY2/2BY2-DD 3BY3/3BY3-DD Integral NaI(Tl) Detector

## Electrical and Mechanical

### Power Required

2BY2 or 3BY3 +12 V at 20 mA. Supplied via a captive power cord terminated in a standard preamplifier power plug (9-pin D connector). Power cord length is nominally 3 m or 10 ft. The preamplifier power plug is compatible with standard preamplifier power connectors provided on most nuclear spectroscopy amplifiers. The preamplifier output signal is delivered on pin-3 of the connector. An optional signal break-out adapter is available for extracting the preamplifier output signal at the power connector.

2BY22-DD or 3BY3-DD Power is provided by the digiDART or digiDART-LF.

### Net Weight

2BY2 or 2BY2-DD 2.5 kg (5.5 lb)

3BY3 or 3BY3-DD 4 kg (8.8 lb)

### Shipping Weight

2BY2 or 2BY2-DD 4.6 kg (10 lb)

3BY3 or 3BY3-DD 4.6 kg (10 lb)

### Dimensions

2BY2 and 2BY2-DD 35.5 cm (14 in.) long x 7.62 cm (3 in.) diameter

3BY3 and 3BY3-DD 38.9 cm (15.3 in.) long x 10.2 cm (4 in.) diameter

## Optional Accessories

**296-ADAPT Signal Break-Out Adapter** Connects to the end of the power cable of the 2BY2 or 3BY3 and separates the preamplifier signal cable from the preamplifier power cable. The adapter 9-pin D connector plugs into the standard preamplifier power connector on the rear of most spectroscopy amplifiers. The 60 cm preamplifier signal cable from the adapter terminates in a male BNC connector for connection to the input of a spectroscopy amplifier.

**C-24-12** 93  $\Omega$ , Coaxial Cable connects the preamplifier output to an amplifier input.

## Ordering Information

| Model     | Description   |
|-----------|---|
| 2BY2      | Integrated 2" x 2" NaI(Tl) Detector   |
| 3BY3      | Integrated 3" x 3" NaI(Tl) Detector   |
| 296-ADAPT | Signal Break-Out Adapter. For use with 2BY2 or 3BY3.  |
| C-24-12   | RG62A/U, 93 $\Omega$ cable with two BNC male plugs; 12-ft. length. For use with 2BY2 or 3BY3. |
| 2BY2-DD   | Integrated 2" x 2" NaI(Tl) Detector for use with digiDART or digiDART-LF.                     |
| 3BY3-DD   | Integrated 3" x 3" NaI(Tl) Detector for use with digiDART or digiDART-LF.                     |

Specifications subject to change  
090309

**ORTEC**<sup>®</sup>

[www.ortec-online.com](http://www.ortec-online.com)

Tel. (865) 482-4411 • Fax (865) 483-0396 • [ortec.info@ametek.com](mailto:ortec.info@ametek.com)  
801 South Illinois Ave., Oak Ridge, TN 37831-0895 U.S.A.  
For International Office Locations, Visit Our Website

**AMETEK**<sup>®</sup>  
ADVANCED MEASUREMENT TECHNOLOGY