


- For counting rate measurements with photons or ion beams
- 10^6 counts/s full scale
- Linear or logarithmic operation
- Fast negative NIM input



Designed for photon or ion beam applications, the ORTEC Model 9349 Log/Lin Ratemeter provides two modes of operation: linear and logarithmic. The linear mode has 11 full-scale ranges from 10 to 10^6 counts/s in 1-3-10 steps. The 5-decade log mode covers the range 10 to 10^6 counts/s in a single span. These choices permit optimum measurement of low, medium, or high steady pulse rates or monitoring of rates that vary through a wide range.

The input signal to the Model 9349 is normally obtained from a discriminator having a NIM-standard, fast-negative logic output signal.

Because of the longer effective scale inherent in its 240° movement and its high accuracy, the unique circular front-panel meter provides excellent readability for both modes.

Zero suppression is provided for up to 100% of any linear range. Any relatively constant background in the counting rate can be subtracted from the data by adjustment of this control. Also, a suppressed zero permits rates that are beyond the nominal full-scale limit to be observed with greater accuracy than could be obtained by switching to a higher range. The choice between 7 linear and 2 log time constants is a further aid to accurate reading of the rate of incoming signals.

In addition to the front-panel meter indications, outputs are provided for both current and voltage recorders, as well as a high-level voltage output for control or monitor applications.

Specifications

PERFORMANCE

LINEAR RANGES 11 ranges from 10 to 10^6 counts/s full scale in 1-3-10 steps.

DEAD TIME <100 ns on the 10^6 range; $<0.3\%$ of average pulse spacing up to the 3×10^4 range; $\leq 1\%$ on the 10^5 and 3×10^5 ranges.

RATED OVERLOAD Maintains full-scale output for X300 overload or 10^7 counts/s, whichever is smaller.

TEMPERATURE INSTABILITY $\leq \pm 0.05\%/^\circ\text{C}$.

NONLINEARITY $\leq \pm 0.15\%$ from 10 to 3×10^4 counts/s range; $\leq \pm 1.5\%$ from 10^5 – 10^6 counts/s.

TIME CONSTANTS 7 selectable time constants, 0.03 to 30 s in 1-3-10 steps.

ZERO SUPPRESSION 0 to 100% of full scale, nonlinearity $\leq \pm 0.25\%$.

LOGARITHMIC RANGE One 5-decade range for 10 to 10^6 counts/s.

TEMPERATURE INSTABILITY $\leq \pm 0.25\%$ of full scale per $^\circ\text{C}$.

ANALOG OUTPUT ERROR $\leq \pm 2.5\%$ of full scale.

STANDARD DEVIATION $\sim 15\%$ with Log Short time constant; $\sim 5\%$ with Log Long time constants.

SLEWING RATE Dependent upon input rate; for any rate change, Log Short time constant provides 10 times faster response than Log Long time constant.

CONTROLS

RANGE 12-position switch selects the full-scale range and either linear or logarithmic mode; linear ranges are 0 – 10 counts/s through 0 – 10^6 counts/s in 1-3-10 steps; log range is 10 – 10^6 counts/s.

TIME CONSTANT 9-position switch selects the time constant for the integrating network; 0.03 to 30 s in a 1-3-10 series for all linear ranges; Short and Long for the log range.

ZERO SUPPRESSION 10-turn precision potentiometer to suppress the zero-reference level for any linear range from 0 to 100% ; the same full-scale span is effective above the preselected zero-reference level.

9349

Log/Lin Ratemeter

INPUT

INPUT Rear-panel BNC connector accepts NIM-standard, fast-negative logic signals, -600 to -1800 mV. $Z_{in} = 50 \Omega$. Minimum pulse width is 4 ns FWHM.

OUTPUTS

PANEL METER 240° circular movement with 8.9 cm (3.5 in.) deflection; accuracy, 2% of full scale; 3 scale markings; 0–1 and 0–3 for linear ranges and 10^{-10}^6 in 5 decades for log range.

ANALOG OUTPUTS BNC connector on rear panel provides 0 to 10 V full scale, dc-coupled with $100\text{-}\Omega$ output impedance.

RECORDER OUTPUTS Binding post connectors on rear panel:

100 mV Provides voltage output with 100 mV full scale; dc-coupled with $100\text{-}\Omega$ output impedance.

1 mA Provides current output of 1 mA full scale; dc-coupled with $10\text{ k}\Omega$ output impedance.

ELECTRICAL AND MECHANICAL

POWER REQUIRED $+24$ V, 50 mA; -24 V, 35 mA; $+12$ V, 30 mA; -12 V, 45 mA.

WEIGHT

Net 1.5 kg (3.5 lb).


Shipping 2.5 kg (5.5 lb).

DIMENSIONS NIM-standard double-width module 6.90 X 22.13 cm (2.70 X 8.714 in.) per DOE/ER-0457T.

Ordering Information

To order, specify:

Model	Description
9349	Log/Lin Ratemeter



Specifications subject to change
011008

ORTEC®

www.ortec-online.com

Tel. (865) 482-4411 • Fax (865) 483-0396 • ortec.info@ametek.com
801 South Illinois Ave., Oak Ridge, TN 37831-0895 U.S.A.
For International Office Locations, Visit Our Website

AMETEK®
ADVANCED MEASUREMENT
TECHNOLOGY