




8-Channel Rapid Analysis Screening of Multiple Gamma Emitting Radionuclides in Clinical Samples



ORTEC Gross Gamma Counter Sample Evaluation Counting System #1				
Sample Date/Time: Friday, October 16 2009 11:13				
	Actual (CPS)	Limit (CPS)	Evaluation	
SampleID: 09-911-000064-1GG Detector: 1	6.36	4.05	FAILED	
SampleID: 09-911-000065-1GG Detector: 2	3.87	3.98	PASSED	
SampleID: 09-911-000066-1GG Detector: 3	3.43	3.80	PASSED	
SampleID: 09-911-000067-1GG Detector: 4	3.54	3.84	PASSED	
SampleID: 09-911-000068-1GG Detector: 5	6.44	3.82	FAILED	
SampleID: 09-911-000069-1GG Detector: 6	3.43	3.68	PASSED	
SampleID: 09-911-000070-1GG Detector: 7				
SampleID: 09-911-000071-1GG Detector: 8				



Easy to Use Rapid and Simple Gamma Radionuclide Screening

- A very compact solution for rapid screening of hundreds of samples/day for possible radioactive contamination.
- No sample processing: Pass/fail evaluations against analysis limits.
- 700 samples may be counted in 20 hrs of operation.
- MDA (LOD)¹ = 52 Bq/L, 300 second count.
- 8 samples (10 ml or 50 ml) counted simultaneously with automated data entry by bar code scanning (optional).
- 8 completely independent NaI(Tl) well detectors and count chains.
- Digitally stable electronics.
- Full spectral data retained for detailed analysis.

¹Minimum Detectable Activity or Limit of Detection.

GammaScreen-8™



8-Channel Rapid Analysis Screening of Multiple Gamma Emitting Radionuclides in Clinical Samples

The post 9-11 world differs drastically from previous nuclear age history. The possibility of a large accidental release of radioactivity affecting a civilian population is now added to by the possibility of an intentional exposure as a result of a terrorist attack involving a radioactive dispersal device or RDD. Following an exposure event or environmental release of radioactive material, large numbers of people may require or demand internal monitoring. Failure to respond quickly can increase the level of public anxiety. Even when the risk is slight, large numbers of subjects may be required to be screened rapidly to reduce public anxiety.

The ORTEC GAMMASCREEN-8 multi-sample screening system has been developed as a simple-to-use tool for this first step screening operation. GAMMASCREEN-8 is a turn-key system based on ORTEC digital signal processing spectroscopy technology.

Eight sodium iodide (NaI (TI)) “well” detectors, each with individual “digiBASE” all-in-one digital spectrometer-photomultiplier bases² are placed within a common and compact lead shield in an arrangement designed to reduce sample cross-talk to a minimum.

Small sample vials may be counted within the well or 50 ml pots may be placed on top of the detectors.




The sample screening software application manages the system ensuring reliable and simple operation. Set-up is very simple; a single parameter entry screen sets the 8 detector bias settings, the stabilization and threshold limits and a common acquisition preset for the sample batches to be screened.

An optional bar code reader can be used to simplify sample ID data entry.

In screening operation, sample gross count-rates are compared to a threshold giving a simple-to-interpret pass-fail criterion. (See example report on previous page).

ORTEC MAESTRO-32[®] MCA software² may be used to scrutinise the detailed spectra of questionable samples; a complete sample load of 8 spectra may be displayed simultaneously, with live spectral display during counting.



“Failed” samples may be investigated further with an HPGe spectrometer, or the sample could be counted longer in-situ for more positive nuclide ID and quantitative analysis using the optional ScintiVision-32™ Quantitative Analysis software².

The ORTEC GAMMASCREEN-8 multi-sample screening system is available as a customized solution to individual requirements.

²More information on the digiBASE digital spectrometer, MAESTRO-32 software, and ScintiVision-32 software is available at www.ortec-online.com.

GammaScreen-8™

8-Channel Rapid Analysis Screening of Multiple Gamma Emitting Radionuclides in Clinical Samples

System Specifications

Lead Shield, floor mounted

Weight 1780 lbs (810kg)

Dimensions (HxWxL) 34.5 x 24 x 28 in
87.6 x 61.0 x 71.1 cm

Width with top open 38.54 in (97.9 cm)

Minimum thickness between detectors 1 in (2.54 cm)

Nal Detectors (8)

Crystal size 3" diameter x 5" long

Sample well 0.67" diameter x 3.9" deep

Sample Tubes (not included but supplied at extra cost if required)

Maximum dimensions 17 mm diameter x 120 mm length with nominal 15 ml volume.

ORTEC digiBase Spectrometer (8)

All-in-one high performance NaI Spectrometer in an integrated photomultiplier base

USB connectivity

Highly stable Digital Signal Processing

Built-in HV supply 0–1200 V

Highly reliable

Low power consumption

100% computer controlled

Minimum Computer Specifications

(individual systems may vary)

Mini-tower format with minimum 1.8 GHz Pentium Dual Core Processor, 1 GB RAM, 160 GB hard drive, 48x32 CD-RW/DVD Combo and 10/100/1000 Gigabit Ethernet.

Windows XP Pro with Vista License.

Office 2007 Pro and Adobe Acrobat 8.1.

USB Hub.

17" LCD Monitor.

HP Laser-Jet 1505 printer or equivalent.

Included Application Software

GS8-B32 GAMMASCREEEN-8 Eight Detector Sample Screening software

A65-B32 MAESTRO-32 MCA Emulation software

Ordering Information

Model

GAMMASCREEEN-8

Description

Complete 8 detector system including all detectors, electronics, PC hardware and software.



GammaScreen-8™

8-Channel Rapid Analysis Screening of
Multiple Gamma Emitting Radionuclides in Clinical Samples

Specifications subject to change
122011

ORTEC®

www.ortec-online.com

Tel. (865) 482-4411 • Fax (865) 483-0396 • ortec.info@ametek.com
801 South Illinois Ave., Oak Ridge, TN 37831-0895 U.S.A.
For International Office Locations, Visit Our Website

AMETEK®
ADVANCED MEASUREMENT TECHNOLOGY