



# High-Purity Germanium (HPGe) Detector Manufacturing

Large single crystals of germanium are grown using the Czochralski technique, which is schematically illustrated by Fig. 4. A precisely cut seed crystal is dipped into the molten germanium and then withdrawn slowly, while maintaining the temperature of the melt just above the freezing point. The rate of crystal withdrawal and temperature of the melt are adjusted to control the growth of the crystal.

Figure 5 shows a crystal during the growth process. High-purity germanium crystals suitable for detector fabrication are almost always grown in a quartz crucible under a hydrogen atmosphere. Near the completion of the growth process, the crystal is tapered gradually at the tail to minimize thermal strain. It is imperative that the crystal be grown to the exhaustion of the melt, because germanium both wets quartz and expands on freezing. The valuable quartz crucible might be fractured if any germanium were left after completion of the crystal growth.

After the crystal is grown and cooled, it is mounted in a Plaster-of-Paris cast for slicing. Figure 6 depicts a mounted crystal during the slicing process. The completed crystal is cut by an ORTEC-designed string saw that causes virtually no damage to the crystal. A slurry of water and silicon carbide is pulled along by a wire, resulting in a sawing action. Sections of the crystal from both top and bottom are checked by Hall effect measurements to determine the impurity concentration and type (n or p). On the basis of the Hall effect results, that part of the crystal which contains detector-grade material is selected. The rejected material is returned to the zone refining operation.

The section of crystal which has both adequate purity and crystallographic perfection for coaxial detector fabrication is then ground perfectly cylindrical. The edge at one end is beveled to a radius ("bulletized") to improve charge collection and timing performance. Figure 7 illustrates the grinding operation. Afterwards, a hole is machined into the unbeveled end so that the central contact of the device may be made later. The detector subsequently is hand lapped all over to remove damage caused by the machining processes.

A lithium diffusion to form the  $n^+$  contact is then performed over the entire outer surface except the flat, unbeveled end for p-type coaxial detectors and on the "walls" of the central hole for n-type coaxial detectors. This lithium-diffused layer is about 600- $\mu\text{m}$  thick. After the lithium diffusion operation, the detector is lapped once more, chemically polished, and a surface protective coating applied. The coating is amorphous germanium hydride deposited by a sputtering process, similar to that described by Hansen, *et al.*, (Ref. 1). Next, the  $p^+$  contact is formed by the ion implantation of boron ions. This last step completes the fabrication process for the coaxial detector element itself. Figure 8 shows schematically the structure of both p-type and n-type coaxial detectors.

At this point the detector is ready to be mounted in a cryostat. The basic function of a cryostat is to cool the germanium detector to its near-liquid-nitrogen operating temperature. For best performance the first stage of the preamplifier is also cooled to low temperature, the entire cold assembly being maintained by the cryostat under high vacuum for both thermal insulation and protection of the internal components from contamination. Figure 9 shows the construction of a typical cryostat system. Cryosorption material (such as selected zeolite or activated charcoal) is used as a residual gas getter or pump to maintain the vacuum for long periods of time. After being loaded into the cryostat, the detector is tested for several parameters, including leakage current and energy resolution. If the device fails a test, it is returned to some previous stage of the process. Another important cryostat design consideration is allowing for the convenient positioning of the detector element. This usually means that a copper cooling rod, needed to conduct the heat from the detector element, is routed from a point outside the dewar to the liquid nitrogen inside the dewar. All of the cryostat materials around the detector should be as low Z as possible to reduce photon scatter. Hence, aluminum, magnesium, beryllium, Teflon, and Mylar are used whenever possible.




Fig. 2. A 3-Coil Zone Refiner.



Fig. 3. A Zone-Refined Ingot.




Fig. 4. Czochralski Technique Schematic.

# High-Purity Germanium (HPGe) Detector Manufacturing

Because many of the steps in the manufacturing process (Fig. 1) have less than a 100% yield, a detector element may spend an extended time in a “loop” before being shipped.

## The Past and the Future

When ORTEC started production of germanium detectors in the late 60’s, a detector of 5% relative efficiency (compared to a 3” x 3” NaI detector) with 3-keV energy resolution at 1.33 MeV was considered typical. In 1975 the “VIP 10” offered 10% efficiency at 2 keV. In 1980 the Ge(Li) technology in use was obsoleted by high purity germanium detectors which could be cycled to room temperature. Recently detectors with 170% efficiency and <2 keV resolution have been produced.

ORTEC pioneered these technologies and others, such as GAMMA-X coaxial detectors, streamline cryostats, and PopTop transplantable detector capsules.

Development work toward larger crystals, application-matched performance parameters, ultra-low background cryostats, and exceptional reliability continues. We welcome suggestions for new, more application-tuned detectors.

## References

1. W.L. Hansen, E.E. Haller, and G.S. Hubbard, “Protective Surface Coatings on Semiconductor Nuclear Radiation Detectors,” *IEEE Trans. on Nucl. Sci.* **NS-27**, No. 1, 1980.




Fig. 5. A Germanium Crystal Being Grown.



Fig. 6. A Mounted Crystal Being Sliced.




Fig. 7. Grinding the Germanium Crystal.




Fig. 8. Schematic of Both P-Type and N-Type Coaxial Detectors.

# High-Purity Germanium (HPGe) Detector Manufacturing




Fig. 9. Exploded View of PopTop Detector Capsule with Horizontal Dipstick Cryostat and 30-Liter Dewar.

Specifications subject to change  
093009

**ORTEC**<sup>®</sup>

[www.ortec-online.com](http://www.ortec-online.com)

Tel. (865) 482-4411 • Fax (865) 483-0396 • [ortec.info@ametek.com](mailto:ortec.info@ametek.com)  
801 South Illinois Ave., Oak Ridge, TN 37831-0895 U.S.A.  
For International Office Locations, Visit Our Website

**AMETEK**<sup>®</sup>  
ADVANCED MEASUREMENT  
TECHNOLOGY