


ORTEC can provide all GEM, GMX, LO-AX, and GWL detectors in low-background configurations, including SGD and PROFILE series GEMS.

The principles of the approach to reducing the background in a HPGe detector operate in two directions: reduction of the background intrinsic to the materials used to make the detector, and secondly reduction of the amount of radiation from any such materials which reach the detector element itself.

### ORTEC Low-Background Detector Materials

- Standard endcaps for low-background GEM and GMX detectors without Be windows are Carbon Fiber. Carbon fiber endcap thickness for low-background detectors is 0.035". Low-background GWL detectors have oxygen-free high conductivity (OFHC) copper endcaps with low-background high purity aluminum well tubes of 0.02" wall thickness.
- If lower energies are required, copper endcaps with low-background Be windows are standard for low-background GMX and LO-AX detectors. The Be window is 0.020" thick for up to 3-1/4" endcaps and 0.030" thick for 3-1/2" endcaps.
- Detector internals (cup, cooling rod, cooling rod clamp and pedestal) are made of OFHC Copper for standard low-background detectors.
- Activated Charcoal is used as the cryo pump, having lower background than sieve.
- Detector flange is made of low cobalt steel.
- XLB option includes a 2 cm thick aged lead back shield between the preamplifier/HV and crystal.



### Low-Background Options for Streamline Detectors

The LB option is constructed as described above.

The XLB option offers an improvement (in most cases) over the LB, by the introduction of the aged lead back shield behind the crystal. (There are instances in which the absence of the lead back shield is preferable, such as when the secondary Pb x rays may create more background at the energy of interest than the additional background present due to the lead back shield's absence.)

### Low-Background Options for PopTop Detectors

The RB option is essentially the same construction as the LB, but the sieve pack in the capsule is retained in order to provide adequate pumping. The detector is shielded from the sieve by an OFHC copper cup.

The RB-B option provides a Be window in a Cu endcap instead of the Carbon Fiber endcap (for GMX and LO-AX detectors only).

### LB Series HJ Cryostat Configuration

The popular HJ configuration is a side-looking cryostat with the preamplifier and high-voltage filter next to the dewar and thus remote from the detector; therefore, the lead back shield immediately behind the detector is unnecessary. Generally, this is the optimum configuration for shielding the detector element from the dewar, the cryogenic pumping material, the preamplifier, and much of the cryostat.

### Factors Affecting Low-Background Gamma Spectroscopy with HPGe Detectors

The radiation background of standard cryostats used by ORTEC for germanium detectors is lower than that required for the majority of users. Net area peak counting rates ~0.1 counts/min are typical at energies of interest. Nonetheless, those measuring environmental samples who require the lowest MDAs in the shortest possible counting time will be best served by a large germanium detector in a LB or XLB cryostat. Coaxial detectors of efficiency from 80% to 175%, with their exceptionally high peak-to-Compton ratios (approaching 100 to 1), are also recommended (Ref. 1).

There are a number of naturally-occurring radionuclides that contribute to the gamma-ray background observed by spectroscopists using a germanium detector system. Contributions from the cosmic-ray induced background, <sup>40</sup>K from building structures, and radon, can be markedly reduced by appropriate shielding (including in some cases an underground location) and flushing the shield with

# Low-Background Germanium Gamma-Ray Detectors

“aged” nitrogen. The principal sources of activity from the cryostat are the primordial emitters,  $^{238}\text{U}$ ,  $^{232}\text{Th}$ ,  $^{235}\text{U}$  and their daughters, man-made radionuclides including  $^{137}\text{Cs}$  fallout, and the activation product  $^{60}\text{Co}$ . There are both full-energy photopeaks and associated Compton background from 46 to 2600 keV.

In some materials, the natural emitter chains may not be in equilibrium. Therefore, ORTEC reports the measured background at all the energies shown in Table 1.

Different spectroscopists have different “low-background” requirements and energies to which the phrase “low-background” applies. Therefore, low-background cannot be rigorously defined as, for example, energy resolution at the 1.33-MeV line of  $^{60}\text{Co}$ . For this reason various laboratories use different nomenclature and different report formats when describing measured background.

To satisfy the needs of spectroscopists, ORTEC does the following:

- Carefully selects low-background materials
- Characterizes completed HPGe detectors in a specialized low-background facility
- Produces the world’s largest detectors in the lowest-background cryostats — the ultimate for gamma spectroscopy of low-level samples
- Provides with each detector an activity report of the 22 most common isotopes
- Remains current with technology for low-background materials. Environmental spectroscopists seeking to minimize MDA while maximizing sample throughput should read: The Benefits of Using Super-Large Germanium Gamma-Ray Detectors for the Quantitative Determination of Environmental Radionuclides.
- For customers that require special configurations, ORTEC can work with the user to supply various materials or utilize customer supplied materials where possible.

## Low-Background Spectroscopy Laboratory and Detector Background Certification Program

Low-background detectors are measured in the ORTEC Low-Background Test Facility, and a certificate is provided reporting the intensity of the lines listed in Table 1.

As a result of the differences in altitude and the extent of radioactivity in materials near to the detector, the background that will be observed at your location may be lower or higher than that measured in the ORTEC Low-Background Test Facility.

The ORTEC Low-Background Test Facility houses two lead shields, configured to accommodate detectors in various configurations. Each shield consists of a 3 inch OFHC copper “well” surrounded by eight inches of low-radiogenic lead. The copper virtually eliminates the lead x-rays resulting from photoelectric interactions of gamma rays with the lead.

For characterizing the materials to be used for detector cryostat construction, the low-background laboratory has a graded-Z (Pb, Cu, Cd) shield containing an ultra-low-background detector. Quality control limits have been established for each type of material dependent on its location in the finished cryostat.

Each completed detector is placed in the appropriate shield and calibrated.

Then, with the source removed, a detector background spectrum is acquired for 100,000 seconds. The background spectrum is searched using a second-difference peak-search algorithm. All identified peaks are visually examined. For any peaks that are part of the list reported as “not found,” a region-of-interest (ROI) is set manually, and the net area is computed by the MCA Emulation software. Finally, a report that lists the ROI net count rates is created.

The data are reported as intensities rather than activities because the activities are a function of the geometry of the calibration, while the intensities are geometry-independent.

The logarithmic plot (Fig. 1), linear plot (Fig. 2), and low-background analysis reports (Tables 2 and 3) show, for comparison, the gamma background spectra of two GEM detectors of identical relative efficiency (56%), one in a standard cryostat and the other in an extra-low-background (XLB) cryostat (Fig. 3). The difference in background is substantial. The background report can be supplied on

Isotope (Parent Nuclide)	Energy in keV	Isotope (Parent Nuclide)	Energy in keV
U x-rays	13.0, 13.3*	$^{137}\text{Cs}$	661.6
$^{231}\text{U}$	25.6*	$^{214}\text{Bi}$ ( $^{238}\text{U}$ )	727.2
$^{137}\text{Cs}$	31.8, 32.2, 36.4*	$^{234\text{mPa}}$ ( $^{238}\text{U}$ )	766.6
$^{210}\text{Pb}$ ( $^{238}\text{U}$ )	46.5	$^{228}\text{Ac}$ ( $^{232}\text{Th}$ )	911.0
$^{234}\text{Th}$ ( $^{238}\text{U}$ )	63.3	$^{228}\text{Ac}$ ( $^{232}\text{Th}$ )	969.0
$^{234}\text{Th}$ ( $^{238}\text{U}$ )	92.6	$^{234\text{mPa}}$ ( $^{238}\text{U}$ )	1001.0
$^{235}\text{U}$ , $^{226}\text{Ra}$	185.7, 186.2	$^{214}\text{Bi}$ ( $^{238}\text{U}$ )	1120.3
$^{212}\text{Pb}$ ( $^{232}\text{Th}$ )	238.6	$^{60}\text{Co}$	1173.0
$^{214}\text{Pb}$ ( $^{238}\text{U}$ )	295.2	$^{214}\text{Bi}$ ( $^{238}\text{U}$ )	1238.0
$^{214}\text{Pb}$ ( $^{238}\text{U}$ )	351.9	$^{60}\text{Co}$	1332.5
Cosmic	511.0	$^{40}\text{K}$	1460.8
$^{208}\text{Tl}$ ( $^{232}\text{Tl}$ )	583.1	$^{214}\text{Bi}$ ( $^{238}\text{U}$ )	1764.5
$^{214}\text{Bi}$ ( $^{238}\text{U}$ )	609.3	$^{208}\text{Tl}$ ( $^{232}\text{Tl}$ )	2614.5

\* The lines lower than 46 keV are reported only for LO-AX and GMX detectors.

# Low-Background Germanium Gamma-Ray Detectors

disk if requested, and can be used in conjunction with the ORTEC MAESTRO MCA Emulation program to visually examine the data or to plot it.



It must be emphasized that the gamma background measured in a detector cannot be better than the background in the laboratory where the detector is operated. For example, health physicists making in vivo measurements should be aware that the beds or chairs in which the subject is placed are generally not made of materials selected for low gamma background and, also, that the subjects being measured emit gamma rays at a rate consistent with the 40–120 nanocuries of  $^{40}\text{K}$  normally found in the human body. Therefore, while every precaution must be taken to obtain good, reliable measurements, it is useless to strive for the “ultimate in low-background” as some physics researchers operating in sophisticated underground laboratories must and can do.

The data obtained at ORTEC are representative of the low-background characteristics of the detector/shield combination in our facility in Oak Ridge, Tennessee. Since background results are dependent upon the shielding, better results may be obtained in sophisticated laboratories (e.g., with active shields or in deep mines).

An example of this for a 40% GEM in a vertical XLB cryostat is shown in Table 4. The third column is from the report of the ORTEC Low-Background Facility. The fourth column contains the reported results on the same detector in a sophisticated low-background laboratory (Ref. 2). Another example is given in Table 5, which reports results obtained at the Gran Sasso National Laboratory, an Italian research facility located under 3500 meters of rock, with a resulting cosmic background reduction of a factor of  $10^6$ . The ORTEC detector has a measured efficiency of 96%. Extraordinary precautions were taken to minimize the gamma background.\*

When comparing the data in Tables 2, 3, and 4, the following considerations must be kept in mind:

As there is no standard, laboratories use different formats for such tables; hence, the obvious differences between Tables 2 and 3 (both obtained at ORTEC) and 4 (obtained at the U.S. National Institute of Standards and Technology).



# Low-Background Germanium Gamma-Ray Detectors

The computer printout reports energy centroids rather than the energy of the nuclides. Therefore, some interpretation is required to understand the centroid/nuclide relationship between Tables 2 and 3 and Table 1. For example, Region #1 in Table 2 (centroid at 45.72 keV) indicates the counts due to <sup>210</sup>Pb in Table 1 (46.5 keV).

When comparing gamma background data obtained from detectors with different efficiency, the difference in efficiency should be factored in, at least in an approximate (linear) way.

A less sophisticated way of characterizing low-background detectors is reporting the total counts per second in a given energy interval, typically from 100 keV to 3 MeV. A large (96% efficiency) ORTEC detector measured at Gran Sasso, the world-class Italian laboratory under a mountain, registered 100 counts per second in that energy interval.

\*C.R. Arpesella, et al., "A Low Background Counting Facility at Laboratori Nazionali del Gran Sasso." (Internal Report LNGS – 92/35 July 1992).

**Table 2. Low-Background Analysis Report for "Standard" Detector.**

```

1      LOW BACKGROUND ANALYSIS
SPECTRUM: A:P40836AN.CHN
ROI FILE: A:P40836AN.ROI
DETECTOR SERIAL NUMBER P40836A; NON LOWBACKGROUND ENDCA  205% RELATIVE EFFICIENCY
ANALYSIS AT 08:13:49 ON 17-AUG-01
MCA NUMBER: 1 SEGMENT NUMBER: 1
REALTIME 100983.20 SECS, LIVETIME 100000.00 SECS
START TIME 13:18:26 ON 09-Jul-99 .0
START CHANNEL  0 LENGTH 8192
SAMPLE MULTIPLIER: 1000.0000
TOTAL CORRECTION FACTOR: 1000.0000
BACKGROUND CALCULATED AS AVERAGE OF  3 CHANNELS
ENERGY = 42.1E-02 + 28.5E-02*CHAN + .0E+00*CHAN**2
    
```

REGION NUMBER	CENTROID ENERGY(keV)	START CHANNEL	STOP CHANNEL	GROSS	NET	PEAK AREA	ERROR NET %	CORRECTED CTS/KSEC
1	46.57	152	172	6442.	643.	25.5	6.43	
2	63.43	210	228	6597.	1065.	14.6	10.65	
3	66.55	230	236	2242.	50.	139.0	.50	
4	76.26	259	281	7245.	219.	84.6	2.19	
5	92.82	313	335	10059.	3036.	6.3	30.36	
6	139.73	477	501	9644.	557.	39.2	5.57	
7	185.94	639	663	10588.	2405.	8.8	24.05	
8	198.59	683	707	8563.	609.	33.5	6.09	
9	238.81	823	849	9801.	1831.	11.7	18.31	
10	277.98	962	988	6393.	318.	57.8	3.18	
11	351.88	1220	1248	5338.	543.	31.1	5.43	
12	510.60	1772	1807	10739.	6233.	3.1	62.33	
13	582.65	2027	2059	4552.	933.	16.8	9.33	
14	608.49	2118	2150	5240.	945.	18.0	9.45	
15	661.15	2302	2332	2990.	0.	.0	.00	
16	668.90	2334	2363	3135.	370.	35.2	3.70	
17	726.38	2532	2566	3153.	219.	65.1	2.19	
1	762.33	2670	2704	3118.	15.	993.1	.15	
19	910.34	3174	3210	2889.	583.	22.5	5.83	
20	961.72	3353	3383	2048.	516.	21.6	5.16	
21	967.65	3385	3411	2045.	330.	29.9	3.30	
22	1000.30	3490	3526	2576.	436.	28.8	4.36	
23	1117.27	3907	3945	2836.	230.	61.3	2.30	
24	1173.92	4091	4129	2055.	86.	142.6	.86	
25	1236.53	4322	4360	2187.	217.	56.3	2.17	
26	1331.28	4650	4686	1384.	0.	0	.00	
27	1459.02	5096	5136	11492.	9659.	1.5	96.59	
28	1762.42	6166	6198	1384.	966.	6.3	9.66	
29	2201.42	7707	7734	688.	282.	18.0	2.82	

END OF ROI LIST

**Table 3. Low-Background Analysis Report for "Low-Background" Detector.**


```

1      LOW BACKGROUND ANALYSIS
SPECTRUM: A:P40836AL.CHN
ROI FILE: A:P40836AL.ROI
DETECTOR SERIAL NUMBER P40836A; LOW BACKGROUND ENDCAP  205% RELATIVE EFFICIENCY
ANALYSIS AT 08:01:57 ON 17-AUG-01
MCA NUMBER: 3 SEGMENT NUMBER: 1
REALTIME 100073.10 SECS, LIVETIME 100000.00 SECS
START TIME 10:24:52 ON 03-Apr-00.1
START CHANNEL  0 LENGTH 16384
SAMPLE MULTIPLIER: 1000.0000
TOTAL CORRECTION FACTOR: 1000.0000
BACKGROUND CALCULATED AS AVERAGE OF  3 CHANNELS
ENERGY = 12.8E-02 + 24.7E-02*CHAN + .0E+00*CHAN**2
    
```

REGION NUMBER	CENTROID ENERGY(keV)	START CHANNEL	STOP CHANNEL	GROSS	NET	PEAK AREA	ERROR NET %	CORRECTED CTS/KSEC
1	45.99	178	194	2045.	116.	75.0	75.0	1.16
2	64.20	250	266	2619.	21.	479.5	.21	.21
3	91.98	363	379	1972.	0.	0	0	.00
4	140.07	559	575	3511.	491.	22.4	4.91	4.91
5	184.31	735	753	3352.	36.	322.4	.36	.36
6	198.56	795	813	3927.	630.	19.0	6.30	6.30
7	236.17	946	964	2710.	0.	0	0	.00
8	278.50	1119	1137	2704.	285.	35.8	2.85	2.85
9	292.57	1173	1191	2123.	17.	546.2	.17	.17
10	348.01	1398	1418	1961.	0.	0	0	.00
11	510.99	2046	2092	8939.	5837.	3.1	58.37	58.37
12	545.70	2324	2346	1244.	2.	3873.7	.02	.02
13	602.81	2429	2451	2045.	0	0	0	.00
14	609.47	2455	2477	2107.	321.	29.5	3.21	3.21
15	657.97	2658	2689	1446.	65.	145.2	.65	.65
16	669.92	2691	2724	1889.	331.	31.3	3.31	3.31
17	727.46	2924	2966	1802.	175.	66.2	1.75	1.75
18	771.64	3084	3126	1714.	23.	518.9	.23	.23
19	910.49	3670	3710	1208.	101.	92.7	1.01	1.01
20	962.29	3883	3911	1351.	534.	13.6	5.34	5.34
21	1001.31	4034	4074	1050.	0.	0	0	.00
22	1120.57	4517	4557	1191.	0.	0	0	.00
23	1173.22	4731	4771	875.	75.	105.4	.75	.75
24	1238.71	4995	5033	856.	108.	69.7	1.08	1.08
25	1332.66	5377	5415	734.	0.	0	0	.00
26	1460.62	5896	5935	1939.	1252.	6.4	12.52	12.52
27	1763.73	7127	7163	455.	159.	30.0	1.59	1.59
28	2613.73	1056	810604	625.	385.	11.9	3.85	3.85

END OF ROI LIST

# Low-Background Germanium Gamma-Ray Detectors



**Table 4. Comparison of ORTEC Low-Background Analysis with that in a Sophisticated Laboratory (NIST).**

LOW-BACKGROUND ANALYSIS  
 DETECTOR SERIAL NUMBER P33P, GEM-40195,  
 Low Background Detector  
 REALTIME 100014.40 secs, LIVETIME 100000.00 secs

CENTROID ENERGY (keV)	NUCLIDE	NET Counts/1000 s (ORTEC LAB.)	NET Counts/1000 s (NIST)
53.6	<sup>73</sup> Ge	1.1	0.4
122.1	<sup>57</sup> Fe	0.33	0.3
145.3	<sup>238</sup> U	0.7	Not detected
186.2	<sup>226</sup> Ra	0.5	Not detected
238.6	<sup>212</sup> Pb	1.3	0.5
295.2	<sup>214</sup> Pb	0.9	0.4
352.0	<sup>214</sup> Pb	1.5	1.0
511.03	β+ annihilation	28	29
583.1	<sup>208</sup> Tl	1.1	0.3
609.3	<sup>214</sup> Bi	1.6	1.0
661.7	<sup>137</sup> Cs	0.6	Not detected
727.1	<sup>228</sup> Ac	0.5	Not detected
1120.0	<sup>214</sup> Bi	0.6	0.3
1173.2	<sup>60</sup> Co	0.3	0.6
1332.5	<sup>60</sup> Co	0.1	0.6
1460.8	<sup>40</sup> K	9	0.3

**Table 5. Counting Rate of Main Gamma Lines for the 96% ORTEC Germanium Detector at Gran Sasso Facility. (1σ errors are indicated.)**

Isotope	Energy	Counts/1000 s
<sup>238</sup> U	295.2	0.22 ± 0.03
	351.9	0.44 ± 0.04
	609.3	0.33 ± 0.03
	1764	0.06 ± 0.01
<sup>232</sup> Th	238.6	0.08 ± 0.02
	583.1	0.04 ± 0.01
	2614.7	0.02 ± 0.01
<sup>40</sup> K	1460.7	0.09 ± 0.02
<sup>137</sup> Cs	661.6	0.05 ± 0.01
<sup>60</sup> Co	1173.2	0.04 ± 0.01
	1332.5	0.03 ± 0.01
<sup>125</sup> Sb	427.9	
	600.6	
	635.9	
<sup>106</sup> Ru	621.8	
	1050.1	

## References

1. R. Keyser, T. Twomey, and S. Wagner, "Benefits of Using Super-Large Germanium Gamma-Ray Detectors for the Quantitative Determination of Environmental Radionuclides," *Radioactivity and Radiochemistry* Vol. 1, No. 2, pp. 46–56 (Spring 1990).
2. R.M. Lindstrom, *et al.*, "A Low-Background Gamma Ray Assay Laboratory for Activation Analysis," *Nucl. Instr. and Meth.*, **A299** (1990) 425–429.

Specifications subject to change  
042910

**ORTEC**<sup>®</sup>

[www.ortec-online.com](http://www.ortec-online.com)

Tel. (865) 482-4411 • Fax (865) 483-0396 • [ortec.info@ametek.com](mailto:ortec.info@ametek.com)  
 801 South Illinois Ave., Oak Ridge, TN 37831-0895 U.S.A.  
 For International Office Locations, Visit Our Website

**AMETEK**<sup>®</sup>  
 ADVANCED MEASUREMENT  
 TECHNOLOGY