

Eichrom Industries provides total solutions for radiochemical separations challenges. Eichrom not only manufactures state of the art extraction chromatography and ion-exchange products, but also provides R&D and consulting services to enhance the value of these products. ORTEC has organized a number of the procedures developed by Eichrom (and its collaborating partners) for convenient access to our Customers, Sales Team, and Technical Support Specialists. For questions about radiochemical separations and alpha-spectrometry please contact:

Michael Schultz, Application Specialist, ORTEC, 801 S. Illinois Avenue, Oak Ridge, TN, 37830-0895, Phone: (865) 481-2446, Fax: (865) 483-0396, Email: michael.schultz@ortec-online.com

## **Technetium-99 in Water (TEVA Disc Method)**

---

### **1. Scope**

- 1.1 This procedure describes a method to separate and measure technetium-99 in water.

### **2. Summary of Method**

- 2.1 Technetium-99 is separated from water samples using Eichrom's TEVA Disc prior to liquid scintillation counting. After concentration of the pertechnetate ion ( $\text{TcO}_4^-$ ) on the disc, Tc-99 is measured by liquid scintillation counting by adding the disc directly to the scintillation cocktail. The sample can be analyzed with and without a Tc-99 spike to determine chemical recovery. Alternately, either of the short-lived gamma emitters, Tc-99m or Tc-95m, can be used as a tracer. The detection limit for this method is 1 pCi/L (37 mBq/L).

### **3. Significance of Use**

- 3.1 This method is a rapid, reliable method for measurement of Tc-99 in environmental samples that combines the selectivity of Eichrom's TEVA Resin with the fast flow rates achievable with a disc format.

#### 4. Interferences

- 4.1 All radionuclides that undergo beta emission or quench the liquid scintillation counting are effectively removed (including C-14, P-32, S-35, Cs-137, and Sr-90) using the Eichrom TEVA disc. Tritium may sometimes follow the technetium due to the absorption of tritium-labeled compounds by the disc. Possible interference by tritium is eliminated by setting the Tc-99 counting window above the maximum energy for tritium beta particles.
- 4.2 Organic matter present in the sample can interfere by quenching during liquid scintillation counting.
- 4.3 Six (6) liters of solution may be passed through an Eichrom TEVA Disc without breakthrough of Tc-99. Larger volumes have not been tested.

#### 5. Apparatus

- 5.1 *Filter apparatus*: suitable to hold 47mm TEVA Disc: Nalgene Reusable Filter Holder recommended (nalgene part # 300-4000 or 300-4050 or 300-4100)
- 5.2 *Forceps*
- 5.3 *Glass beakers*
- 5.4 *Liquid scintillation counter*
- 5.5 *Liquid scintillation vials, glass*
- 5.6 *Watch glass*

#### 6. Reagents

- 6.1 Unless otherwise indicated, all references to water should be understood to mean deionized distilled water.
- 6.2 *Hydrogen peroxide (30 wt %)*
- 6.3 *Liquid scintillation cocktail* - Ultima Gold<sup>®</sup>-LLT (Packard Instrument Company)
- 6.4 *Nitric acid solution (0.01M)*- Add 640  $\mu$ L concentrated HNO<sub>3</sub> (sp gr 1.42) to 900 mL water and dilute to 1L with water.

6.5 *TEVA Disc* - 47 mm, Part number: TE-D50-FX (50 discs) or TE-D200-FX (200 discs)

## 7. Procedure

### 7.1 *Water Sample Preparation:*

- 7.1.1 Measure the sample volume using a graduated cylinder (or equivalent) and transfer the volume to an appropriate size beaker.
- 7.1.2 Adjust the pH of the sample to pH 2, if required.
- 7.1.3 Analyze the sample with and without adding Tc-99 spike to determine chemical recovery.

Note: An alternative is to use Tc-99m or Tc-95m as a tracer, measuring the gamma activity of Tc-99m or Tc-95m. If Tc-99m is used, allow the Tc-99m to decay for approximately 1 week after measuring the gamma activity and before measuring Tc-99 using liquid scintillation counting.

- 7.1.4 Add 10 mL of 30 wt % H<sub>2</sub>O<sub>2</sub> (per liter of sample) to the beaker, stir and cover with a watch glass. Remove cover once the sample begins to boil.
- 7.1.5 Heat the sample to about 90°C for 1 hour to oxidize Tc(IV) to Tc(VII), forming TcO<sub>4</sub><sup>-</sup>, oxidize some of the organics present, and destroy excess H<sub>2</sub>O<sub>2</sub>.
- 7.1.6 If bubbling due to decomposition of the hydrogen peroxide has not stopped as the sample cools, continue heating until bubbling has stopped. Stir occasionally with glass stirring rod.
- 7.1.7 Allow the beaker to cool to room temperature.
- 7.1.8 If the sample contains insoluble matter, filter the sample to remove solids prior to using the TEVA Disc.

### 7.2. *Eichrom TEVA Disc separation*

- 7.2.1 Using forceps, carefully place a TEVA Disc on the filtering apparatus.

- 7.2.2 Place the funnel over the disc and transfer the water sample from step 7.1.7 into the filtering funnel.
- 7.2.3 Allow the sample solution to filter through the disc by gravity or vacuum. (No vacuum is necessary. Depending on which filtering apparatus is used, it is possible to achieve a flow rate as high as 33 mL/min with gravity. An average of 98% Tc-99 recovery achieved using 200 mL/min with vacuum on 1 L simulated ground water sample type)
- 7.2.4 After filtering the sample, rinse the original beaker or container with the minimal volume of water required (depending on beaker size) and transfer this rinse into the filtering funnel.
- 7.2.5 Allow the rinse solution to drain through the disc completely.
- 7.2.6 Rinse disc with 50 mL 0.01M nitric acid solution.
- 7.2.7 Squirt 5 mL of water using a squeeze bottle to rinse the sides of the funnel. This will remove any residual material of the disc which may have suspended during filtration.
- 7.2.8 Connect the filtering apparatus to a vacuum pump ( if gravity flow used to this point) and apply vacuum for a few seconds to remove any residual liquid from the disc. Alternately, if the filtration apparatus has a Luer fitting, attach a 50 cc syringe and draw residual liquid from the disc.

### 7.3 *Counting preparation:*

- 7.3.1 With forceps, remove the TEVA Disc from the filtering apparatus and transfer it into a liquid scintillation vial.
  - Note: The disc can be rolled gently while it is still on the filter apparatus so that it can be easily inserted into the LSC vial.
- 7.3.2 Add 15 mL of Ultima Gold LLT into the LSC vial containing the disc. Cap the vial and vortex until the disc disintegrates and forms a colorless translucent gel.
- 7.3.3 Set up the scintillation counting window to measure from 26 Kev to 350 Kev or according to lab protocol.

## 8. Calculations

### 8.1 Calculate chemical yield

#### 8.1.1 Measurement of spiked and unspiked samples

$$Y = \frac{(C_S - C_U)}{E \times A_S}$$

where:

- Y = chemical yield
- C<sub>S</sub> = count rate spiked sample, cpm
- C<sub>U</sub> = count rate unspiked sample, cpm
- E<sub>S</sub> = counting efficiency for Tc-99
- A<sub>S</sub> = Tc-99 tracer activity, dpm, corrected for decay from reference date

#### 8.1.2 Tc-99m (or Tc-95m) Tracer Method

$$Y = \frac{(C_S - B_S)}{E_S \times A_S}$$

where:

- C<sub>S</sub> = measured Tc-99m (or Tc-95m) tracer, gamma cpm
- B<sub>S</sub> = background, gamma cpm
- E<sub>S</sub> = gamma counting efficiency for Tc-99m (or Tc-95m)
- A<sub>S</sub> = Tc-99m (or Tc-95m) tracer activity, dpm, corrected for decay from reference date

### 8.2 Calculate the Tc-99 activity as follows:

$$A = \frac{(S - B)}{E \times V \times Y \times 2.22}$$

where:

- A = Tc-99 activity in the sample (pCi/L)
- S = sample counts/time in minutes, cpm

- B = blank counts/time in minutes, cpm  
E = counting efficiency for Tc-99 determined from a direct spike of Tc-99 in LSC vial containing a TEVA Disc dissolved in 15 mL Ultima Gold LLT.  
V = sample volume, L  
Y = chemical yield (determined in section 8.1)

### **REFERENCE**

- (1) Eichrom Analytical Procedure, TCW01, rev. 1.3, "Technetium-99 in Water," March 1, 1995.
- (2) Sullivan, T., et al. "Determination of Technetium-99 in Borehole Waters using an Extraction Chromatographic Resin." 37th Annual Conference on Bioassay, Analytical and Environmental Radiochemistry. Ottawa, Canada. 1991.
- (3) Mann, Darrin. Interview: K-25 Plant, Oak Ridge Site.
- (4) DOE Methods Compendium. Report # RP550. "Technetium-99 Analysis Using Extraction Chromatography."
- (5) Wyse, E.J. and Fadeff, S.K., "Alternative Techniques for the Determination of Technetium-99 in Groundwaters: ICP/MS and Extraction Resin." Submitted for publication.