

1.0 Introduction

MAESTRO Version 7 introduced a standardized method for collecting spectrum data using the List Mode capability provided in ORTEC MCAs (specifically the digiBASE, digiBASE-E, DSPEC Pro, DSPEC 50/502, and IDM-200). The List Mode files (.LIS) can also be loaded into a MAESTRO buffer window and filtered by acquisition time to isolated spectrum data to a specific time range. This capability is further described in the MAESTRO literature and Software User's Manual.

Although standard functionality for reading, filtering data by date range, and saving List Mode files is fully supported in MAESTRO, some users may want to read the List Mode files with other applications. This document describes the format of the List Mode files specifically for this purpose.

Notes Regarding Real and Live Time when using List Mode:

- 1) The file content of the .LIS file can be slightly different depending on how the acquisition is terminated. When data acquisition is stopped using the STOP command (menus, toolbar, or Job File) the Read/Live times are read immediately and no further List data is read. The file header data is then updated, and the file is closed. The instrument, however, may take a short time (up to a second or so depending on the connection speed and instrument type) before the acquisition actually stops. Since MAESTRO continues to update the sidebar with the Real and Live time until the instrument indicates that the acquisition has completed, there will often be a small difference between the Real and Live time displayed on the detector window display as compared to the Real and Live time stored in the .LIS file. The Real and Live times stored in the .LIS file correspond with the actual list data from the instrument.
- 2) When collecting List Mode data with the digiBASE the Live Time is always set identical to the Real time when .LIS files are recalled. This is the result of the Live Time information not being stored in the instrument when operating in this mode. If dead time is critical for the measurement, then data collection in PHA mode or a different instrument choice should be considered.

2.0 File Structure

The List Mode file consists of a 256 byte header followed by a series of 32-bit List Mode data records. The List Mode data record format is dictated by the information available in each instrument as described in the sections below.

2.1 List Mode Header Record

Offset	Size (Bytes)	Data Type	Description
0-3	4	32-bit integer	List Mode header format (must be -13)
4-7	4	32-bit integer	List data style 1=digiBASE, 2=PRO List, 3=Not Used, 4=digiBASE-E
8-15	8	double	Acquisition start time in OLE DATE format
16-95	80	char(80)	Device address as a string (suitable for CONN control address)
96-104	9	char(9)	MCB Type string (SHOW_VER response)
105-120	16	char(16)	Device serial number as a string
121-200	80	char(80)	Text description of data
201	1	char(1)	Non-zero if energy calibration constants are valid
202-205	4	char(4)	Null terminated text string representation of energy units (e.g. "keV")
206-209	4	float	Offset (lowest order) coefficient in the energy calibration polynomial
210-213	4	float	Linear (middle) coefficient in the energy calibration polynomial
214-217	4	float	Quadratic (highest order) coefficient in the energy calibration polynomial
218	1	char(1)	Non-zero if the shape calibration constants are valid
219-222	4	float	Offset (lowest order) coefficient in the shape calibration polynomial
223-226	4	float	Linear (middle) coefficient in the shape calibration polynomial
227-230	4	float	Quadratic (highest order) coefficient in the shape calibration polynomial
231-234	4	32-bit integer	Optional conversion gain (Valid if non-zero)
235-238	4	32-bit integer	Optional Detector ID number (valid if non-zero)
239-242	4	float	Optional Real Time in seconds (Valid if non-zero)
243-246	4	float	Optional Live Time in seconds (Valid if non-zero)
247-255	9	char(9)	Unused (Round to 256 bytes for the header)

2.2 List Mode Raw Data Format

The List Mode data consists of a series of 32-bit data records that are formatted based on the List Data Style specified in the second parameter of the Header Record. The following sections describe the data record formats for each List Data Style

2.2.1 digiBASE Format (List Data Style = 1)

The digiBASE transmits one 32-bit word for every detectable event in the NAI detector as well as other events. The 32-bit words are encoded as follows:

Bit	Description
31	0 for event
30-21	Amplitude of pulse
20-0	Time in microseconds that the event occurred

The time is a 21-bit number in units of microseconds. The number rolls over to 0 every 2.097152 seconds. In order to track the rollovers, a “time only” event is sent from the digiBASE to the computer every 1.048576 seconds. The format of the “time only” event is as follows:

Bit	Description
31	1 for time-only
30-0	Current Time in microseconds

Note that Live Time is not recorded by the digiBASE. When a digiBASE format file is read by MAESTRO, the Live Time is set to the Real Time and dead time shows on the sidebar as zero.

ORTEC MAESTRO List Mode File Formats

2.2.2 PRO List Format (List Data Style = 2)

This format is applicable to the DSPEC Pro, IDM-200, and DSPEC 50/502 instruments.

These instruments transmit one 32-bit word for every detectable event in the detector as well as other events. The 32-bit words are encoded as follows:

	31	30	29		16	15		0	
ADC Word	1	1	14-bit ADC value			16-bit real time in 200 nS ticks			
RT Word	1	0	30-bit real time in 10 mS ticks						
LT Word	0	1	30-bit live time in 10 mS ticks						

	31		24	23		16	15		8	7	0
Hdw Time	00000000			00000000			16-bit real time in 200 nS ticks				
UMCBI Time	00000001			Byte 2			Byte 1		Byte 0		
UMCBI Time	00000010			Byte 5			Byte 4		Byte 3		
UMCBI Time	00000011			00000000			Byte 7		Byte 6		
ADC CRM	00000100			00000000			Counts per 10 mS period				
Ext Counter 1	00000101			00000000			Counts per 10 mS period				
Ext Counter 2	00000110			00000000			Counts per 10 mS period				
GM Counter (IDM-200 Only)	00000111			00000000			Counts per 10 mS period				

The ADC Word contains a raw 14-bit ADC value and a 16-bit real time stamp. The Hardware Time Word contains the same 16-bit real time stamp. In both cases this value represents the time that the word was created as represented by an internal clock that increments every 200 nS. The 200 nS clock rolls over from 49999 to 0, which gives it a period of 10 mS.

Every time the 200 nS clock rolls over a RT Word and LT Word are created. The 30-bit real time stamp on the RT word corresponds to the number of times the 200 nS clock rolls over. The 10 mS real time clock and the 10 mS live time clocks are cleared when the CLEAR command is sent to the spectrometer. Otherwise they increment every 10 mS while data is collecting. Using the RT Word as a time standard and the 200 nS time stamps from the ADC word you can determine the time that any ADC event occurred, ± 200 nS, relative to the beginning of the acquisition.

LT Words are created by the spectrometer at the same time as are RT Words. The 30-bit live time stamp represents the ADC live seconds since the beginning of the acquisition (or since the last CLEAR command). By calculating the relative change in live time between

ORTEC MAESTRO List Mode File Formats

any two LT Words in a List Mode stream you can calculate the live time for a given block of data, ± 10 mS.

The ADC CRM word contains the value of the spectrometer's count-rate meter and is generated at the same time as the RT and LT words (i.e. it always appears next to these words in the List Mode data stream). The count-rate meter is incremented each time the ADC fast channel detects a pulse. Since many pulses may be rejected before one is converted into an ADC event, this number is not the simple sum of the ADC events over the past 10 mS but instead it represents the true input count-rate, as well as it can be measured by the instrument. Note that this value is zeroed after an ADC CRM word is created therefore it represents only 10 mS of ADC input counts.

The Ext Counter words contain the value of the external input pulse counters and are generated at the same time as the RT and LT words. The external pulse counters count the positive pulses at the external input of the spectrometer. Note that the external inputs must stay high for at least 40 nS for the pulse to be counted. These counters are zeroed after an Ext Counter word is created therefore it represents only 10 mS of external pulses.

Optionally, the external counters may be configured to be a simple level sensitive input. In this case the counter value is either 0 or 1 depending on the level of the input at the time the counter word is generated.

2.2.3 digiBASE-E Format

The digiBASE-E transmits one 32-bit word for every detectable event in the detector as well as other events. The 32-bit words are encoded as follows:

ADC Word	1 1	17-bit time pre-scale in 80 nS Ticks	13-bit ADC channel number
RT Word	1 0	30-bit Real Time counter in 10 mS Ticks	
LT Word	0 1	30 Bit Live Time counter in 10 mS Ticks	
Ext Sync	0 0	13-bit Real Time counter in 10 mS Ticks	17-bit time pre-scale in 80 nS ticks

The Real Time timestamp (in mS) is calculated by multiplying the last real time record by 10 and the Live Time timestamp (in mS) is calculated by multiplying the last live time record by 10.

The ADC event time is calculated (in mS) by adding the last Real Time (in mS) value to the pre-scale value multiplied by $8E-5$. (where $8E-5 = 80 \text{ nS ticks} * 1E-6 \text{ mS/nS}$)

The Ext Sync words contain the value of the external input pulse counters. The external pulse counters count the positive pulses at the external input of the spectrometer. The sync time is calculated by adding the real time counter to the time pre-scale value.