

Reports the Way You Need Them from a Searchable Database

- Custom Report Generation from GammaVision[®]-32
- Multiple custom report formats via specific templates
- Spectral plot included with report
- Creates Access[®]-compatible results database
- Selective searching of results database
- Archival capability

Reports Your Way Report Writer (A44-B32) in conjunction with GammaVision-32 lets you produce a wide range of output reports — in exactly the form you require.

Instantly Searchable Database Report Writer reads in the analysis results from the GammaVision-32 output file, stores them into an Access-compatible database, and produces the report. Thereby, you progressively generate a PC-searchable database of all your analysis results. Then use Access to quickly generate, for example, daily summaries, exception reports, and cumulative release reports. Report Writer thus acts as a bridge between GammaVision-32 results and database applications or a laboratory information management system.

Our Template or Yours The output report is generated using a template file. The default template file supplied reproduces the standard GammaVision-32 output report. Additional user-defined template files¹ may be used to produce totally custom output reports via this template development tool.

In addition to automatic output to the printer, the reports may be saved to disk in various industry standard formats: ASCII, Rich Text, Crystal Reports, WORD, Comma-separated.

Multicolor Plots Automatically A plotting program is included which permits the spectrum to be plotted to the printer automatically as part of the output report. The plot attributes may be prespecified.

As Automatic as You Want to Be Report Writer can be run independently in order to regenerate or reformat reports from previous analyses, or automatically as part of the GammaVision-32 analysis sequence, in which case the standard internal GammaVision-32 output report is disabled.

Good Housekeeping A facility is provided to allow database and supporting files to be archived to a specified directory. File types may be chosen from templates, spectral files, output files, and region-of-interest files. The archival process requires only a single click. Similarly a “purge” option allows removal of such files to the system recycle bin or permanent deletion.

Ordering Information

Model Number	Description
A44-B32	Report Writer option for GammaVision-32
A44-U32	Update to A44-B32
A44-G32	Documentation for A44-B32
A44-CR	SEAGATE Crystal Reports Professional

Prerequisites

A44-B32 will function satisfactorily on any system capable of supporting GammaVision-32 .

For custom report template development, a copy of Crystal Reports Version 7.0 or later is required.

¹ In Crystal Report Writer, V7.0 format; trademark of Crystal Software, Vancouver, BC.