

Model	DSP	NIM	Stabilizer	Interface	Conversion Time	Channels	Special Features
ASPEC-927	No	Yes	Yes	USB 2	<2.0 μ s	16k	Dual MCA, Zero Dead Time Correction
926-M32-USB	No	Yes	No	USB, Parallel Printer Port	8 μ s	8k	
920E*	No	Yes	No	Ethernet or USB**	15 μ s	4k	16 input (shared ADC)
921E*	No	Yes	Yes	Ethernet or USB**	1.5 μ s	16k	High count rate for analog systems
TRUMP-PCI-2K	No	No	Yes	PCI card	8 μ s	2k	
TRUMP-PCI-8K	No	No	Yes	PCI card	8 μ s	8k	

*Requires purchase of MAESTRO-32 (A65-B32) Software.

**Requires purchase of the model DPM-USB and software drivers.

Specifications subject to change
022111