

These charts are intended as selection guides only. For complete and precise specifications, consult the data sheet for each model.

Model	449	661	9349
Type	Log/Linear	Linear	Log/Linear
Full Scale Ranges Minimum Maximum	10 counts/s 10 ⁶ counts/s	25 counts/s 10 ⁷ counts/s	10 counts/s 10 ⁶ counts/s
Input Types	±Logic Pulses, >3 V high and >50 ns wide	Positive discriminator; Negative NIM	Negative NIM >4 ns wide
Pulse-Pair Resolution	<100 ns, or <1% of average pulse spacing	<40 ns	<100 ns, or <1% of average pulse spacing
Audible Output	449-2 option	No	No
Recorder Output	Yes	Yes	Yes
Package and Width	NIM-2	NIM-1	NIM-2
Special Feature	0 to 100% zero suppression	Fast Response Circuit	0 to 100% zero suppression

Specifications subject to change
082609